

Smart Protein Layers **Real-time Normalization of Western Blot Analysis**

Jana Heise^{1,2}, Jan Heise¹


¹NH DyeAGNOSTICS GmbH, Weinbergweg, 23 06120 Halle, Germany ² Weinberg Laboratories, Weinbergweg 23, 06120 Halle, Germany

The Problem of Western Blot Normalization

Western blot analysis is a method to detect a target protein by immunochemistry in a complex sample. The technique includes protein separation, its transfer to a membrane and diverse washing and antibody incubation steps.

Reliable comparison of detected target signal intensities between different samples requires an appropriate way of normalization.


Reported fold changes of the target protein must not be based on an artifact of a reference signal (e.g. saturated pixel) and/ or experimental errors.



Normalization

Lane	1	2	3	4
Signal intensity px (AU)	<u>.</u>	<u>.</u>	<u>.</u>	<u>.</u>
Target	1.890	2.230	3.480	3.610
Reference 1 (Loading control 1)	2.350	12.670	5.320	4.810
Reference 2 (Loading control 2)	< 65.000	< 65.000	< 65.000	<65.000
normalized Target / Ref. 1	1,00	0,22 ?	0,81 ?	0,93 ?
normalized Target / Ref. 2	1,00	n.d.	n.d.	n.d.

The Solution



additional information about the sample, enable automated data evaluation.


		And a state of the
		Construction of the local division of the lo
		Concession of the local division of the loca
		And in case of the local division of the loc
		Contraction of the local division of the loc
		And Address of the Owner, where the Owner,
		and the second se
a strengthere and		
		100 C
and the second second		1000
		And Inc. of Concession, name
and the second second		tion - Side
		1000
and the second second		THE OWNER WHEN THE OWNER
and the second s		Construction of the


Signal intensity px (AU)	,			
Target	1.890	2.230	3.480	3.610
Reference real-time total	42.640	44.240	35.770	37.650
normalized Target	1,00	1,14	1,54	1,67

Conclusions

SPL offers

- real-time total protein detection in gels and blots
- normalization of the protein content between samples
- step-by-step monitoring of the complete Western Blot workflow
- precise normalization of the target protein expression in Western Blots
- accurate comparison of target protein expression between different experiments





SP real-time gel & blot monitoring

