Smart Protein Layers

Quick Guide Protein Labeling

1 Required materials for SPL Labeling

- Smart Label working solution (SLW = tube marked with „Smart Label Reagent B)
- SPL Buffer
- SPL Smartalyzer (SMA)
- 60 mM DTT (freshly prepared)
- Calibrator Mix
- Protein sample (Protein concentration: max. 10 µg/µl)

2 Preparation of the SLW (only before first usage)

1. Allow vials containing Smart Label Reagent A and B to warm up to room temperature and spin down briefly.
2. To solubilise the Smart Label add 15 µl of Smart Label Reagent B (contains molecular sieves to avoid water pollution) to Smart Label Reagent A and mix. Spin down briefly.
3. Transfer all liquid from Smart Label Reagent A to Smart Label Reagent B, mix and spin down briefly.
4. The SLW solution (=“Smart Label Reagent B”) is now ready to use and stable for 6 month. Store the SLW solution at -20 to -80 °C.

3 Preparation of the Calibrator (CAL) as a Master Mix is recommended

1. Transfer 8 µl CAL A per gel into a fresh micro-centrifugation tube.
2. Add 2 µl CAL B* per gel and mix (a dilution might be necessary to use with low abundant targets).
3. Add 2 µl of 60 mM DTT per gel (freshly prepared), mix and centrifuge briefly.
4. Denature the proteins by heating the mixture for 5 min at 95 °C.
5. Apply 12 µl of the prepared CAL to each gel.

*CAL B is not needed for usage of SPL Gel Kits.

4 Preparation of the Reaction and Loading Mix (RL Mix)

6 µl	SPL Buffer	}	x (n**+1)	**no. of samples
2 µl	60 mM DTT			
2 µl	SMA (S or L)			

Note:

SMA as well as CAL B can be used in dilution when the signal is too prominent. Please pre-dilute SMA und CAL B with water or PBS before preparing the RL or CAL Mix (e.g. 1:100). The volumes used for RL and CAL Mix should not be changed.

5 SPL Labeling

Variant A) Equal volumes with different protein concentrations (only possible with SPL)

protein conc.	sample vol.	H ₂ O***	RL Mix	SLW	denaturing	gel load
up to 10 µg/µl	10 µl	-	10 µl	1 µl	5 min, 95 °C	21 µl

Variant B) Equal Protein concentrations using different volumes (e.g.)

50 µg	4 µl	6 µl	10 µl	1 µl	5 min, 95 °C	21 µl
50 µg	9 µl	1 µl	10 µl	1 µl	5 min, 95 °C	21 µl

***or buffer in which sample is solved

For an excel sheet ask the NHD service team (info@dyeagnostics.com).

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