




The VELUM Dry Blotter

Easy as pie protein transfer in HD quality, buffer-free

The VELUM Dry Blotter allows for a very simple high quality, high efficiency protein transfer from gels to PVDF and nitrocellulose membranes. Compared to existing blotting devices the handling is much easier as its special combination of pressure and capillary action does not require any buffers or electric power. To ensure the homogeneity of the applied pressure the VELUM Dry Blotter comes with a torque key. Blotting area: 265 x 120 mm.



The VELUM Dry Blotter



VELUM Dry Blotter
Product. No. PR88
List price: 895,00 EUR

Specialprice*
599,00 EUR

For orders or further
information please contact

NH DyeAGNOSTICS GmbH
Fon: +49 (0)345 2799 6413
Fax: +49 (0)345 2799 6412
Email: info@dyeagnostics.com

*Excluding tax, price EXW (Incoterm 2000)

www.dyeagnostics.com